


- * **Immaculate in performance**
- * **Immaculate in appearance**
- * **Immaculate in construction**



* The LEEC Precision Incubator operating between 0 C. and 90 C. has the most reliable heating stability to be found in any modern incubator. At 37 C. the temperature variation is less than 0.5 C. throughout the entire cabinet.

* Here it is, the modern styled Incubator. Immaculate in design and appearance.

* The construction of this Incubator has been designed and it is built of materials to give many years of unflinching service. All the component parts are readily accessible.


















Write to the makers for details of the LEEC range of Incubators, and ask for information about the LEEC self-contained Refrigerator Unit which operates with the LEEC Incubator (with cooling coil) at temperatures from ambient to 0 C.

LEE**C**

makers of quality laboratory apparatus

LABORATORY & ELECTRICAL ENGINEERING CO.

71-73 Goldsmith Street Nottingham Tel: Nottm. 48331

| | | |
|--|--|---|
|  <p>"E.R.L." HOLDERS for TEST TUBES, AMPOULES, BOTTLES, etc. In 8 sizes either in black synthetic rubber, or softer red organic rubber</p> |  <p>"ESCORUBBER" RUBBER and SILICONE RUBBER BUNGS In various qualities and many sizes</p> |  <p>"ESCORUBBER" ADJUSTABLE HOLDER for TEST TUBES & WATER BATH RACK Various sizes of holes available at 20/- each</p> |
|  <p>"ESCOPLASTIC" VINYL TUBING Non-pyrogenic and non-toxic for Blood Transfusion, sterilizable at 120°C., flexible, excellent for general laboratory and hospital work, milk, various chemicals, hot water, etc.</p> | |  <p>"ESCORUBBER" SILICONE RUBBER TUBING Non-toxic, translucent, will withstand temperatures up to 200/250°C.</p> |
|  <p>"ESCORUBBER" RUBBER and SILICONE RUBBER BOTTLE CAP LINERS Made in many standard sizes or to specification</p> |  <p>"E.R.L." RUBBER FILTERING CONES for BUCHNER FLASKS In sets of 6 at 22s. 6d. per set</p> |  <p>"ESCORUBBER" RUBBER and SILICONE RUBBER SHEETS Made in various qualities and many sizes</p> |
|  <p>"E.R.L." SILICONE RUBBER COVERED MAGNETS for MAGNETIC STIRRING 1" at 7s. 6d. each 1½" at 9s. 6d. each 2" at 12s. 6d. each</p> |  <p>"ESCORUBBER" SR 300 SILICONE RUBBER MOULDING PASTE</p> |  <p>"E.R.L." RUBBER PIPETTE TEATS Extra heavy, in 7 colours for identification purposes approx. 2½ c.c. and 5 c.c. also approx. 7½ c.c. & 10 c.c. (red only)</p> |
|  <p>"ESCORUBBER" R 33 RUBBER TUBING High quality for various laboratory and hospital applications</p> |  <p>"ESCORUBBER" CLOSURES for BLOOD and SALINE SOLUTION BOTTLES Ref. No. P2—made in high quality rubber, and used extensively by many well-known hospitals</p> |  <p>"ESCOPLASTIC" PLASTIC PIPETTE TEATS With remarkable resistance to many chemicals In 7 colours for identification purposes at 45s. per gross</p> |
|  <p>"ESCOPLASTIC" POLYTHENE DISPOSABLE GLOVES 1000 at 16s. 8d. per 100 plus 1s. 4d. P.T.</p> |  <p>"ESCOPLASTIC" PLASTIC DISPOSABLE PETRI DISHES</p> |  <p>STAINLESS STEEL GRADUATED MEASURES and OTHER LABORATORY and HOSPITAL EQUIPMENT</p> |
| <p style="text-align: center;">ESCO (RUBBER) LTD. Manufacturers of Rubber, Silicone Rubber and Plastic Products SEAL WORKS · SEAL STREET · LONDON · E.8</p> | | |

WHAT'S NEW FROM B.D.H.?

ALGINIC ACID 40-100 MESH

(B.P.C. 1963, Appendix 8)

The special grade of alginic acid referred to in the B.P.C.¹ for separating certain B group vitamins by the ion-exchange chromatographic method has been prepared and tested in the B.D.H. laboratories, using Foster and Murfin's² modification of Specker and Hartkamp's method³.

1. British Pharmaceutical Codex 1963, Appendix 8, p. 1291
2. Foster, J. S. and Murfin, J. W., *J. Pharm. Pharmacol.*, 1961, **13**, Suppl. 126T-132T
3. Specker, H. and Hartkamp, H., *Z. anal. Chem.*, 1953, **140**, 167-70
10 g 33/-

N-ACETYL-L-CYSTEINE *A mucus digestant for microbiological examinations*

A reagent for digesting and homogenising sputum samples, of particular value in diagnosis and for chemotherapeutic investigations into tuberculosis treatment. The most tenacious sputum is liquefied in 5-30 seconds¹, permitting the use of a more dilute sodium hydroxide decontaminant, and giving appreciable increases in the number of acid-fast organisms surviving the initial treatment². This novel application of *N*-acetyl-L-cysteine is described in a leaflet, sent free on request.

1. Kubica, G. P., *et al.*, *Amer. Rev. Resp. Dis.*, 1963, **87**, 775-9
2. Webb, W. R., *J. Thor. Cardio. Surg.*, 1962, **44**, 330
25 g 22/6

SODIUM CHLORIDE SOLUTION 18% w/v *for preparing isotonic (0.9% w/v) solution*

A concentrated saline solution from which 0.9% w/v isotonic solution for laboratory use may be prepared by simple 20:1 dilution with water; a convenient and time-saving alternative for workers who usually make up their own solution by weighing and dissolving sodium chloride.

1 litre 8/6 2½ litres 17/3

'TEEPOL' 610 *for serum iron determinations*

When 'Teepol' 610, an ionic detergent based on sodium sec.-alkyl sulphates, is used in determinations of serum iron, the initial precipitation of proteins is unnecessary. Colorimetric methods employing this technique have been reported by Webster¹, who estimated serum iron after intravenous injections of iron-dextran, and, more recently, by Sandford². The method can be applied to both fresh and frozen sera, and to the determination of chelated iron in urine. ('Teepol' 610 is manufactured by the Shell Chemical Co. Ltd., and distributed for laboratory purposes in the United Kingdom by The British Drug Houses Ltd. 'Teepol' is a trade mark of the Shell Chemical Co.).

1. Webster, D., *J. Clin. Pathol.*, 1960, **13**, 246-8
2. Sandford, R., *J. Clin. Pathol.*, 1963, **16**, 174-7
1 litre 9/6



for proven accuracy, precision and reliability —



SGOT, SGPT

DADE TRANSAMINASE REAGENTS are used in a simple, direct, and proven procedure for colorimetrically measuring the levels of serum glutamic oxalacetic transaminase (SGOT) and serum glutamic pyruvic transaminase (SGPT). Each set contains stable, balanced reagents for 100 determinations.

DADE TRANSAMINASE REAGENTS are standardized to yield the same values as obtained by the U. V. reference method of Karmen and Wroblewski. You prepare only one calibration curve and verify it with Enza-trol®, the known enzyme test control serum.

LDH REAGENTS

DADE LACTIC DEHYDROGENASE (LDH) REAGENTS are used in the determination of serum or body fluid LDH levels. Elevations of LDH are often demonstrated in cases of infarctions or malignancies. The test is conveniently performed in a vial containing DPNH.

LDH and SGOT tests, when run in parallel, provide valuable additional information as an aid in the diagnosis of myocardial infarctions

HBD REAGENTS

DADE HYDROXYBUTYRATE DEHYDROGENASE (HBD) REAGENTS, for the quantitative colorimetric measurement of HBD activity, have a diagnostic advantage over determinations of total serum LDH and SGOT in cases of myocardial infarction. The peak activity of HBD is obtained at the same time as LDH but remains elevated for a longer period than either LDH or SGOT. HBD has been reported to be specific for indicating heart tissue damage.

HBD determinations, when performed in conjunction with the LDH test, can differentiate between liver damage and heart tissue damage.

ENZA-TROL®

DADE ENZA-TROL® fulfills the need for a known enzyme test control serum. It provides you with a reliable check on substrates, reagents, instrument variation and technic. Used exactly as serum, Enza-trol® contains serum protein and known values of enzymes normally found in human serum. Because most enzymes are unstable in solution, Enza-trol® is supplied in the dried form. After reconstitution, stability is the same as enzymes in human serum.

DADE

AUTHORISED DISTRIBUTORS IN THE UNITED KINGDOM
ARMOUR PHARMACEUTICAL COMPANY LTD / Hampden Park / EASTBOURNE / Sussex / Eastbourne 51111

ELECTRICALS (LABORATORY) LTD / 'Hedgefield House' / Blaydon-on-Tyne / Co DURHAM / Ryton 2381

ELECTRICALS (LABORATORY) LTD / 3A St Vincent Street / EDINBURGH 3 / CAL 2188

ARNOLD R HORWELL LTD / 17 Cricklewood Broadway / LONDON NW2 / GLA 8777 & 8631

KERNICK & SON LTD / Monthermer Road / CARDIFF / South Wales / Cardiff 26128



ORME SCIENTIFIC LTD / Bentley Avenue / Middleton / MANCHESTER / Middleton 2648

REPUBLIC OF IRELAND

E. J. BROWNE, M.P.S.I.
 6 The Crescent / Sallinoggin / Dunloaghaire / DUBLIN / EIRE / Dublin 807020

— you can depend on Dade

Reagents for
coagulation tests



**Deprothrombinised
Ox Plasma**
(as used for OWRENS
P&P method)

Price: £8 per 100 vials
each containing
2 ml Ox Plasma
£4 per 50 vials each
containing 2 ml
Ox Plasma
£2 per 25 vials each
containing 2 ml
Ox Plasma

**Rabbit Brain
Thromboplastin**
(Acetone Dried)

Stability proven.
Results Reproducible.
Price: £1 15 0d. per box
of six vials—each vial
sufficient for 5ml
of Thromboplastin

**Rabbit Brain
Thromboplastin**
(Saline Extract—
Phenolized)

Price: £1 15 0d per box
of 6x5 ml bottles.

**Bell & Alton
Platelet Substitute**
(For Thromboplastin
Generation Tests)

Price: 18/- per box
of 6x5 ml bottles.




BRITISH PRODUCTS
OF BRITISH RESEARCH

Stayne Laboratories Ltd.,
Hillbottom Road, Sands,
High Wycombe, Bucks.
Tel. High Wycombe 5721

When "FREEZING" is the only way

Certain techniques (for example: enzyme techniques, fluorescent antibody techniques, or where rapid diagnosis is required) prohibit the use of fixatives during the preparation of tissue sections for microscopy. Because freezing keeps sections fresh, unfixed sections of frozen tissue are in increasing demand. Although special freezing cabinets can be constructed to enclose the microtome, the CAMBRIDGE MICROTOME CRADLE has the great advantage that it can be used with certain commercially available deep-freeze cabinets already installed in many laboratories (without modifications). The bottom of the cabinet is still available for normal storage when the microtome is in position. OPERATION OF MICROTOME:—Once in position, the microtome can be operated directly, without complicated mechanisms, through openings (rubber sealed) at the top. The effective design of this cradle has been confirmed by two years' routine operation of prototype models.



Cambridge

MICROTOME CRADLE for deep freeze operation

MICROTOMES—THE CAMBRIDGE RANGE

- FOR ROUTINE OPTICAL MICROSCOPY *The Cambridge Rocking Microtome*
- FOR LARGE SECTIONS *The Cambridge Rotary-Rocking Microtome*
- FOR ROUTINE ELECTRON MICROSCOPY *The Cambridge Fine-Section Microtome*
- FOR ADVANCED WORK—SPECIALIST APPLICATION *The Cambridge Ultra-Microtome (A. F. Huxley pattern)*



Write for further information, quoting CAB/AY. For future publications ask for Mailing form AY/7/64

CAMBRIDGE INSTRUMENT COMPANY LTD. 13 GROSVENOR PLACE, LONDON, S.W.1

KOCH-LIGHT LABORATORIES LTD.

Have pleasure in announcing the publication of:

"ENZYMES, COENZYMES and SUBSTRATES 1964/5"

A completely revised and extended Catalogue giving a comprehensive listing of Bio-chemicals with purities, activities and numerous literature references

"NEW COMPOUNDS"

An interim Catalogue listing over 600 new products

"CUSTOM SYNTHESIS"

A leaflet describing our custom synthesis facilities

"METABOLIC PATHWAYS CHART" —3rd edition

A revised edition of this well-known chart

"PENTEX BLOOD PROTEINS"

A completely new listing of these Pentex products

Write or phone Dept. P.I.

KOCH-LIGHT LABORATORIES LTD.

Telephone: COLNBROOK 2262

BUCKINGHAMSHIRE