



# **Detect Rheumatoid Factor and Thyroid Autoantibodies by rapid slide tests**



**Rheuma-Wellcotest\***

(Code No KT 05)

**Thyro-Wellcotest\***

(Code No KT 06)

**Kits**

Each kit utilizes a sensitized polystyrene latex suspension and contains all the necessary materials and full details for performing the slide tests and interpreting the results. The Rheuma-Wellcotest may also be used for quantitative estimations.



**Wellcome Reagents Limited**  
Wellcome Research Laboratories  
Beckenham Kent England BR3 3BS

Trade Mark\*



THE STANDARD PROCEDURE FOR BOTH QUALITATIVE AND QUANTITATIVE  
ESTIMATIONS OF HUMAN CHORIONIC GONADOTROPHIN IS

THE **PREGNOSTICON** TEST

## **NOW ORGANON SIMPLIFY THE STANDARD PROCEDURE**

### **PREGNOSTICON 'ALL-IN'**

is based on the same haemagglutination principles  
as the Pregnosticon test but it has been modified so that all  
the reagents are freeze-dried in one snap-top ampoule which also  
acts as the reaction tube.

1. "A rational development towards the  
simplification of existing methods"  
(D., Munch. med. Wschr. 111 (2) 100-  
102, 1969).

2. A simpler technique — one ampoule  
contains all test materials.

3. Eradication of preparatory procedures.

4. "Economy in materials and apparatus  
was among the technical advantages  
noted" (J. Med. Lab. Technol. 27, 492-  
494, 1970).

## **PREGNOSTICON 'ALL-IN'**

### **From ORGANON**



Further information available on request.  
Organon Laboratories Ltd · Crown House · Morden · Surrey

---