



Thrombo-Wellcotest*

Latex slide test for detection of fibrin (ogen) degradation products in serum and urine.

Really complete and at the ready!

Bottles stand upright for ease of identification and accessibility.

20 test kit : ordering Code number HA 13.

Sensitivity : is controlled at 2µg/ml and carefully balanced to minimise the possibility of false positives.

Specificity : latex particles coated with antibody raised against fragments D & E. Therefore, sensitive to both early and late FDP.

Note the full count

- 3.0 ml latex suspension
- 1 ml positive control serum and
- 1 ml negative control serum
- 2 x 25 ml glycine saline buffer
- 20 vacuum draw FDP collection tubes
- 1 tube plastic pipettes
- 1 tube mixing sticks
- Reaction tile (6 rings)
- 1 rubber teat
- Direction circular

Kit shelf life : 12 months (2°-8°C)

**The kit with the confidence . . .
... to stand up and be counted**

**Thrombo-Wellcotest—
Supreme in FDP detection**



Wellcome Diagnostics
A Division of The Wellcome Foundation Limited Temple Hill Dartford England DA1 5AH

Please supply me with further details on Thrombo-Wellcotest

Name _____

Position _____

Address _____

Wellcome Diagnostics
A Division of The Wellcome Foundation Limited,
Temple Hill, Dartford, England, DA1 5AH.

* Trade Mark

