

Treatment with Kabikinase within hours of the onset of serious occlusion usually prevents irreversible tissue necrosis and may often be limb-saving. Kabikinase therapy should be commenced as soon as possible using the simplified standard dosage scheme.\* No elaborate assays are necessary, but as thrombotic damage usually proceeds rapidly after occlusion it is advisable to hold in the Pharmacy standby stock for at least one patient treatment. Otherwise delivery is arranged by the swiftest transport available at the time from Kabi Pharmaceuticals Ltd. Telephone 01-567 4717. \*Brit. med. J., 1966, i, 454

(pure streptokinase for intravascular use)
IN THROMBOSIS AND EMBOLISM
Manufactured by AB Kabi, Stockholm, Sweden

Full details, reprints and literature from:



## **BRIGHT'S** — THE INTERNATIONAL NAME FOR CRYOSTATS ANNOUNCE 7 OUTSTANDING NEW MODELS

THE ONLY CRYOSTATS APPROVED BY THE PATHOLOGICAL DEPT. OF THE ROYAL COLLEGE OF SURGEONS OF ENGLAND

## BRIGHT'S STAR FFATURES



BENCH TYPE MODELS BM/CS BM/FCS BM/RS

- New self-aligning Anti-Roll device fitted with micro adjustment.
- Quick Action Knife and Object Holder Freezer.
- Microtome Chamber Stainless Steel.
- Rollers titted for easy movement of cabinet.
- Window free from frost and mist
- Freezing Chamber free from cooling coil obstruction therefore easy to
- Control Panel in full view of operator.
- Internal fluorescent light completely recessed.
- All Models designed to be operated from a comfortable sitting position.




FREE-STANDING MODELS FS/CS



BENCH MODEL designed for use where floor space is unavailable to accommodate the Tree-Standing

accommodate

CRYOSTATS fitted with Rotary Rocker or Modified Rocking Cambridge Microtome Operated by remote controls. For every type of research or routine microsectioning



manufactured from rustproof steel and finished with grey Stets its, permanent PVC surface, resists scratches, acids, detergents and time BRIGHT'S MODEL ES TESIS III. with a Thermo Electric cooled micro-

ALL BRIGHT'S CRYOSTATS as

tome installed in a conventionally cooled environment. Knile, tissue and environment temperatures all adjustable

BRIGHT'S NEW OT US CRYO STAT is compactly designed and cabe operated from a comfortable sitting position. Ideal in operating theatres and autopsy rooms. Fee-nomic in price, space and connect

Full particulars sent on application

## REFRIGERATOR SERVICE LTD. **BRIGHT'S**

Cryostat Division, Factory No. 3, Clifton Road, Huntingdon, Hunts. Telephone: Huntingdon 4528 Telegrams: Cryomatic Huntingdon