


NEW VERIFY CITRATE VERIFY I VERIFY II


extended-time control plasmas for coagulation tests



Verify is a control plasma containing reduced amounts of those clotting factors which are depressed by anticoagulant therapy: Factors II, VII, IX, and X. During anticoagulant therapy the prothrombin activity is maintained between 10% & 30% of Normal. **Verify I** approximates to plasma with 30% prothrombin activity, and **Verify II** approximates the 10% level.

Verify is more sensitive than diluted normal plasma to minute variations in test conditions. It is also simpler to use – only one pipette operation is necessary for reconstitution. Verify is versatile, and may be used as an 'abnormal control' in detecting bleeding disorders as well as in monitoring anticoagulant therapy.

Full information will be provided on request

Albumin Kabi	human
Ceruloplasmin Kabi	human
Fibrinogen Kabi	human & bovine
Gamma Globulin Kabi	human
Haptoglobin Kabi	human
Plasmin Kabi	human
Plasminogen Kabi	human
Transferrin Kabi	human
Manufactured by AB KABI Stockholm Sweden	
 <p><i>For information and prices write to sole distributors in U.K.:</i> Kabi Pharmaceuticals Ltd. Bilton House, Uxbridge Road, Ealing, London W.5. Telephone 01-567 4717</p>	