

**Simplified procedure**

**NOW**

Results in 70 minutes from primary plate  
or four hour broth culture

**NOW**

Rock slide, read result in one minute or less

**NOW**

Definitive result from  
fewer organisms

**NOW**

37°C incubation

**NOW**

No centrifugation

**Streptex**

makes the decision in  
Streptococcal grouping  
still with

**Confidence  
Coverage  
Certainty**

Excellent correlation with Lancefield  
Comprehensive kit, covering A, B, C, D, F and G  
Highly-specific latex reagents giving clear-cut results



**Wellcome Reagents Limited**  
299-303 Hither Green Lane  
London SE13 6TL, England

\* Trade Mark

JCP 1

Please send me your detailed Streptex brochure.

Name \_\_\_\_\_

Institution \_\_\_\_\_

Address \_\_\_\_\_



**Wellcome Reagents Limited**  
299-303 Hither Green Lane  
London SE13 6TL, England


# Just give us half a minute



In just 30 seconds the Coulter® S7 Automated Haematology Analyser provides the standard 7 parameter profile from whole or pre-diluted blood.

Using the proven Coulter® principle, the compact Model S7 counts the red and white cells simultaneously in two separate chambers for greater accuracy and precision. It incorporates a new electronically controlled diluter, easier calibration with digital readout, and can dispense 10ml of diluent to simplify capillary blood dilutions, using only 44.7µl of specimen.

The S7 is designed to be fast, reliable, accurate, easy to use and is economical to purchase and maintain.



## Coulter Electronics

Coulter Electronics Limited,  
Northwell Drive,  
Luton,  
Bedfordshire, LU3 3RH  
England.  
Telephone: (0582) 582442  
Telex: 825074

**COULTER® - count on us**