

ID sample	Abbott Real-Time PCR		Idylla SARS-CoV-2 Test						RealStar SARS-CoV-2 Real-Time PCR		
	Ct value	Result	Cq value (N)		Cq value (Orf1b)			Result	Ct value (E)	Ct value (S)	Result
1	30,54	Positive	39,18	38,41	41,06	39,51	41,54	Positive	31,46	31,67	Positive
16	24,83	Positive	38,28	38,13	38,37	39,81	39,29	Positive	28,12	30,45	Positive
25	35,00	Positive					17,44	Positive	25,78	26,31	Positive
30	26,00	Positive	35,10	35,56	37,31	36,65	37,32	Positive	29,84	27,78	Positive
31	12,00	Positive	31,06	30,98	32,67	32,7	32,71	Positive	14,24	22,69	Positive

Supplementary Table 2. High discrepant results between Idylla SARS-CoV-2 and Abbott Real-Time PCR were further analyzed by a third technique (RealStar SARS-CoV-2 Real-time PCR).