

ID sample	Abbott Real-Time PCR		Idylla SARS-CoV-2 Test					RealStar SARS-CoV-2 Real-Time PCR			
	Ct value	Result	Cq value (N)		Cq value (Orf1b)		Result	Ct value (E)	Ct value (S)	Result	
1	37,70	Uncertain	43,02	42,68	40,97		40,96	Positive	35,18		Uncertain
2	38,95	Uncertain	42,47					Negative			Negative
3	38,86	Uncertain						Negative			Negative
4	37,95	Uncertain	41,15	41,36		40,17		Positive	35,64		Uncertain
5	37,37	Uncertain	41,26	41,82	41,43			Positive	33,54	34,06	Positive
6	38,88	Uncertain						Negative			Negative
7	38,76	Uncertain	42,25	42,66			43,40	Positive	36,18	33,78	Positive
8	39,38	Uncertain		42,15				Negative			Negative
9	38,98	Uncertain						Negative			Negative
10	39,64	Uncertain						Negative			Negative
11	37,94	Uncertain		42,67	41,45			Positive			Negative
12	38,28	Uncertain		43,45				Negative			Negative
13	37,83	Uncertain						Negative			Negative
14	37,90	Uncertain	41,84	42,07				Negative			Negative

Supplementary Table 3 Comparison among Idylla SARS-CoV-2, Abbott Real-Time PCR and RealStar SARS-CoV-2 Real-Time PCR in Abbott Real-Time PCR uncertain cases (Ct>37).